

Handling alcohol-based surface disinfectants



Pre-soaked Tissues

For immediate use.



Ready-to-use Solution*

Apply liquid with a disposable wipe.



Ready-to-use Foam

Spray foam onto the surface or onto a disposable wipe.

Application

CORRECT:

- ✓ Remove visible dirt from surfaces before disinfection
- ✓ Wet surfaces evenly with liquid
- ✓ Wipe with evenly exerted pressure
- ✓ Wait until surfaces have dried completely
- ✓ Switch off ignition sources and ventilate well, as products with high alcohol content are highly flammable
- ✓ After use, dispose of wipes in household waste
- ✓ Close packaging immediately after use

WRONG:

- ✗ Dilute disinfectant
- ✗ Disinfect floors or larger surfaces with alcohol-based products
- ✗ Use for hand disinfection
- ✗ Damp wiping is enough
- ✗ Wipe the surfaces dry
- ✗ Allow disinfectant to get inside the device

Products

for rapid disinfection against bacteria, yeasts and enveloped viruses e.g. HIV, influenza viruses and coronaviruses



Products with high alcohol content

Only use on alcohol-resistant surfaces (e.g. not on acrylic or plexiglass)



Pre-soaked Tissues:

- Bacillol® Tissues
- Bacillol® AF Tissues

Ready-to-use Solution

- Bacillol® AF

Products with low alcohol content

Can also be used on sensitive plastics, acrylic and plexiglass as well as monitors, displays, keyboards, etc.



Pre-soaked Tissues:

- Bacillol® 30 Tissues

Ready-to-use Foam:

- Bacillol® 30 Foam

* German jurisdiction: max. 50 ml per m² of surface in accordance with TRGS 525

Use surface disinfectants safely. Always read the label and product information before use.



Research for
infection protection
bode-science-center.com

Download
Poster:

